



MEDISEPT

Perfect disinfection



MEDISEPT Wipes for hands and surfaces disinfection

High-quality wipes for disinfecting hands and all small alcohol-resistant surfaces

Efficient wipes for disinfecting hands and surfaces, including surfaces of medical devices. It has bactericidal, yeasticidal, mycobactericidal, tuberculocidal, and limited virucidal effect (e.g. Herpes Simplex, HIV, Hepatitis B and C, Rota virus). Dermatologically tested.

Product features:

- Tested in accordance with new standards for disinfectants.
- Hygienic disinfection in 30 sec. (soaking liquid tested in acc. with EN 1500).
- High quality, non-woven fabric: 50 g/m².
- The capacity of one wipe is 0.5 m with a double rubbing.*.
- Wipe dimensions: 19 x 15 cm.

* The test was carried out based on the methodology of the EN 16615 standard for a four-field test, on a porous surface, at 24 ° C

User benefits:

- Recommended for hands and small areas in the medical area and in public places.
- Convenient in usage.
- Delicate, fresh fragrance.
- Contain a moisturizer - glycerine.
- The closing window made of hard plastic guarantees a constant moisture after opening.



Active against
Tbc, Rota
in 30 sec.



50 pcs packaging



Disinfect
and protects
hands



Disinfects
surfaces

MEDISEPT Wipes for hand and surfaces disinfection

Application:

MEDISEPT Wipes are high-quality alcoholic wipes designed for hygienic hands disinfection and any small alcohol-resistant surfaces. They are recommended for surfaces such as tabletops, knobs, handles, seating surfaces, medical equipment, and accessories made of plastic, stainless steel, porcelain and ABS plastic. Due to its intended use on the skin, the product contains care-substance additives which moisturize the hands.

Effectiveness	EN14885	Time
Bacteria (including MRSA)	EN 13727	30 sec.
Bacteria, fungi (C. albicans)	EN 16615	30 sec.
Fungi (C. albicans)	EN 13624	15 sec.
M. Tuberculosis (M. terrae)	EN 14348	30 sec.
Enveloped viruses (HIV, HBV, HCV, Vaccina, Influenza, Vaccinia, BVDV, Herpes simplex, Ebola, SARS-Cov-2)	EN 14476	30 sec.
Rota Virus	EN 14476	45 sec.
Hygienic disinfection*	EN 1500	30 sec.

* test performed on saturated formulation,
1 tissue contains approx. 3.56 ml of disinfecting liquid

Directions for use:

Take a wipe out of the package and wipe hands or surface. Close the package tightly. After opening, use the wipes within 21 days. Do not use if the wipe is dry. Do not use on areas sensitive to alcohol.

Ingredients: 50 pieces of tissues contain:

178 g of saturated formulation

100 g of formulation contains:

60 g - Propan-2-ol,
0.5 g - d-gluconic acid, compound with N,N''-bis(4-chlorophenyl)-3,12-diimino-2,4,11,13- tetraazatetradecanediamidine (2:1)

Packaging unit:



50 wipes in a flow pack with self-adhesive closing

Permission number for biocidal products in Poland: 6678/16

Medical device class II within the meaning of the Medical Devices Act.
Intended for disinfecting non-invasive medical devices and for hygienic hand disinfection.

The safety data sheet and the component data sheet are available at www.medisept.pl

CE 2274

 Producer
MEDISEPT Sp. z o.o.
Konopnica 159 c
21-030 Motycz

phone: + 48 81 535 22 22
export@medisept.pl
www.medisept.pl

Revision: 04.2020