

For more information or orders, please fill the form below and send it back to the address mentioned on the back side.

- Name:

- First name:

- Street:

- Zip code:

- City:

- Country:

- Telephon No.:

- Doctor:

- Comment:

.....

.....

----- Information

----- Order ----- VPMR-System
(please mark)

Signature:

Reference-No.:

VPMR-System	Manual-powered vacuumpump with accessories
EM-250	lubrication jelly
EM-B 100	Constriction-band set with ringloader and lubrication jelly
EM-M 5	Constriction-band, size 5
EM-M 6	Constriction-band, size 6
EM-M 7	Constriction-band, size 7
EM-M 8	Constriction-band, size 8

Contact:

EuroMedical GmbH
Ligusterweg 2
83313 Siegsdorf
Germany
Tel: +49 8662 7050
Fax: +49 8662 668549
E -mail: info@euromedical-gmbh.de
www.euromedical-gmbh.de



presented by:

VPMR SYSTEM
Vacuum therapy



with
erectile
dysfunction

Simple and quick help with impotence.
 Sure, successful, affordable ... for all ages.

VACUUM THERAPY is the easiest way to solve the problem of impotence without side effects.

VACUUM THERAPY is safe, easy and reasonable in price. Doctors prescribe this therapy successfully since more than 25 years. Contra-indications are not known. This device is used by millions of men in Europe and USA.

VACUUM THERAPY functions by a mechanically generated vacuum, achieving a good erection. A successful unification with your partner is thus possible, without having to resort to any surgical intervention, injections or tablets.



The vacuum therapy system with EuroMedical® can be used without restrictions.

A plastic cylinder, the penile tube, is placed over the penis, upto the penis root and up against the body. Vacuum is generated and blood flows into the penis making it grow and stiffen, an erection which makes a sexual unification possible. An elastic ring, the constriction/ tension band, is placed over the penis to prevent the venous blood from 'leaking back' from the penis. Once the penis is stiff and the constriction band is properly positioned, the cylinder can be removed.



Impotence can result from diseases like Diabetes mellitus or Arteriosclerosis, appears after surgery, spine or internal injuries and after use of various medications. Impotence can be hormone-caused or result from smoking, alcoholism or drug abuse.

Benefits of vacuum therapy are

- a quick erection
- a significant strengthening of virility
- short-term commencement
- use without risks
- a determinable time of erection from up to 30 minutes



VPMR System

consists of

- a unit with manual vacuum pump and penile cylinder
- several constriction/tension bands and adaptor bushings
- a loading cone
- lubrication jelly
- carrying case