

DISINFECTANTS

AIRDECON[®] 200

ADVANCED WHOLE ROOM DECONTAMINATION SYSTEM (WRDS)

The airdecon[®] 200 is a Whole Room Decontamination System (WRDS), providing professionals with a highly effective means of Infection Prevention & Control. The airdecon 200tm delivers a dry mist of Ozone and Hydrogen Peroxide for decontaminating air and surfaces. It is fitted with a unique diffusion principle that allows the disinfectant to be in contact, by displacement, with surfaces being treated, including those which appear inaccessible.

The device is fitted with a complex dynamic dispersal system, including a fluid vaporising system and an embedded ozone generator, all together linked to an HCI (Human Computer Interface). Used in conjunction with Amity HP 75 (H2O2) or peroxone (O3+H2O2). It is a fully automated system.

For patient wards, emergency rooms, operating rooms and infectious disease units.

For room decontamination, just roll it in place. Simple and quick, keeps your facility working.

Primary features attributable to airdecon[®] 200 are:

- Simple and quick - keeps your facility working
- Safe & effective - Log7 reduction in one operation
- Peroxone - unique dosing system for high performance
- Ozone generation and detection
- Hydrogen Peroxide dry mist
- Amity HP 75 stored in sealed and ready-to-use containers with RFID
- Non-corrosive
- Bio-degradable
- Complies with NFT 72-281 (2009)



In the Box

- Airdecon[®] 200
- An IEC standard power cord
- LAN (RJ45) cable
- CD-Rom (traceability program included)
- WIFI & tablet computer available as cost option

DISINFECTANTS

Application

The airdecon® 200 was developed to be used either Amity HP 75 [H₂O₂] or peroxone (ozone + Amity HP 75) [O₃+H₂O₂].


The peroxone process involves two essential steps: ozone dissolution and hydrogen peroxide addition.

There are two major effects from the coupling of ozone with hydrogen peroxide:


Oxidation efficiency is increased by conversion of ozone molecules to hydroxyl radicals; and

Ozone transfer from the gas phase to the liquid is improved due to an increase in ozone reaction rates.

Antimicrobial data confirmed to EN 13704, EN 14347, EN 14348, EN 14563, EN 13727, EN 1040, EN 14561, EN 14476, EN 13624, EN 14562 & NFT 72-281 (2009).



Decontaminaton



Directives HARMONIZED STANDARDS
2006/42/EC NF EN 61010-1
2002/96/EC Electrical and Electronic Equipment Directive
2002/95/EC ROHS Directive
2004/108/EC EN 55022
CEI 61000-4-2
CEI 61000-4-3
CEI 61000-4-4
CEI 61000-4-5
CEI 61000-4-6
CEI 61000-4-11
CEI 61000-3-2
CEI 61000-3-3

AIRDECON® 200

WHOLE ROOM DECONTAMINATION SYSTEM (WRDS)

airdecon® 200 properties

- 300m³ capacity.
- System ensures the device does not start or automatically stops running if movement is detected and triggers an alarm.
- Ozone generation and detection.
- Humidity sensor.
- 7.5% Hydrogen Peroxide (Amity HP 75).
- RFID modules ensure that only authorized/ recognized containers can be inserted.
- Calibrated temperature sensor triggers an alarm when temperature exceeds 50°C.
- Clear warning lights show operators when decontamination is in preparation.
- 5 to 8 µm particle size.
- Filtering of all disinfectant to avoid the risk of nozzle blockage and disinfectant contamination.
- No pump, no liquid pressure sensor installed.
- Amity HP 75 stored in sealed and ready-to-use containers.
- No risk of leakage or accidental spillage (safe for operator).



Amity Limited, Libra House, West Street,
Worsbrough Dale, Barnsley, South Yorkshire,
S70 5PG, UK

T: +44 (0) 1226 770787

F: +44 (0) 1226 770757

E: sales@amityinternational.com

WWW.AIHC.CO.UK

© Amity Limited all rights reserved 2018

DISTRIBUTED BY: