

Compact HbA1c Analyzer

PocketChem A1c Advanced



Reliable

Fast & accurate result CV3%, in 4.5min.

User friendly

LCD touch screen
Multi-lingual voice & video guide
Built-in Thermal Printer

Easy to use

Single test cartridge containing synthetic materials (Bead, Washing solution)







PocketChem A1c is a fully automated analyzer utilizing boronate affinity method for determination of HbA1c in human whole blood.


Simple test steps







HbA1c Test cartridge



- Easy blood sampling
- No sample preparation
- •No pipetting required

Related items





Test strip for urine albumin and creatinine AUTION SCREEN Microalbumin / Creatinine



For early detection of chronic kidney disease

PocketChem A1c Test Kit

Within-Run Precision (IFCC sample)

Each of the samples was measured 20 times Per each for precision using 2 meters across 3 lots

IFCC calibrator	Normal	Abnormal
Mean HbA1c %	5.5	8.9
CV (%)	2.4	1.9

Day-to-day Precision (IFCC sample)

Each of the samples was measured 1 time a day for 20days using 2 meters across 3 lots.

IFCC calibrator	Normal	Abnormal
Mean HbA1c %	5.4	9.0
CV (%)	2.6	1.9

Specification	
Measured sample	Capillary whole blood/venous
Measurement principle	Boronate affinity
Measurement range	4.0~15.0%(20~140mmol/mol)
Required sample volume	3.5uL
Measurement time	4.5min
Operating condition	Temperature: 17-32℃
	Humidity: 15~75% R.H.

^{*}Design and specification may be changed without prior notice.

Printer type	Internal printer
Shelf life	12 months(Reagent)10 cartridge/ box
Out put	USB port
Size	205w×200d×150h
Weight	1.6Kg
Option	External Barcode Reader(option)
Storage condition	2~30°C(Reagent)

arkray global business, inc.

Yousuien-nai, 59 Gansuin-cho, Kamigyo-ku, Kyoto 602-0008, JAPAN

TEL 81-75-662-8979 FAX +81-75-431-1202