


ALLERG-O-LIQ

Specific IgE REAST for the quantitative determination of specific IgE in human serum or plasma with biotinylated allergens

The ALLERG-O-LIQ is based on the “reversed enzyme allerge sorbent test (REAST)” and represents an innovative, flexible and highly accurate test system using anti-human IgE coated microtiter plates in combination with biotinylated reagents (Figure 1). Compared to conventional IgE detection systems all interfering serum components are removed from the system by washing. This leads to the advantage of just measuring allergen-specific IgE and finally to the high sensitivity and specificity of the test.

Allergen-specific IgE is the mediator of immediate hypersensitivity (Type I allergy), which is often accompanied by increased titers of total IgE. Therefore, the quantitative in-vitro determination of specific IgE and total IgE represents an important diagnostic tool for patients with suspected Type I allergy.

spec. IgE REAST REF 0520960FL	Σ	96
spec. IgE REAST REF 0524800FL	Σ	480



slgE REAST Specifications

- ▲ „Reversed“ ELISA technology
- ▲ Several hundred single allergens and allergen mixes incl. recombinant and nativ, highly-purified allergens are available
- ▲ Determination of specific and total IgE in the same test run
- ▲ Manual and fully automated test procedure possible
- ▲ Incubation time: 3 hours at RT
- ▲ Detection system: HRP/TMB ($OD_{450nm/620nm}$)
- ▲ 6 point calibration according to WHO Standard
- ▲ Wide measuring range (0,35-100 IU/mL) according to RAST classes 0-6 for slgE and 5-1000 IU/mL for total IgE

Addition of stop solution and reading of the optical density at 450 nm
 * Incubation time of the substrate for manual procedure is 30 min

Figure 1
Principle of the ALLERG-O-LIQ

Performance against IVD methodes

Agreement between ALLERG-O-LIQ and ImmunoCAP® (ThermoScientific) was investigated using 495 specific IgE values to Dermatophagoides pteronyssinus (d1), Dermatophagoides farinae (d2), timothy grass pollen (g6), birch pollen (t3) and hazel pollen (t4), cat dander (e1) and mugwort pollen (w6). Receiver Operating Characteristic (ROC) analysis and spearman correlations were performed for all seven allergens together. An excellent area under the curve (AUC) value was found at 0.97 (Confidence interval, CI 0.96-0.98) for all seven different allergens (see figure 2) compared to ImmunoCAP® results.




Figure 2
ROC analysis for ALLERG-O-LIQ vs. ImmunoCAP® for seven different allergens (d1, d2, g6, t3, t4, e1 and w6) with n=495 results.

Spearman correlation between ALLERG-O-LIQ and ImmunoCAP® for all seven allergens reveals a correlation coefficient of 0.91 (CI 0.90-0.93, see figure 3).




Figure 3
Spearman correlation between ALLERG-O-LIQ and ImmunoCAP® for seven different allergens (d1, d2, g6, t3, t4, e1 and w6) with n=495 results.

Performance against skin prick test and nasal Provocation

Skin Prick test (SPT) results combined with nasal provocation (NP) results were compared to the ALLERG-O-LIQ. AUC value of ALLERG-O-LIQ compared to 287 SPT and NP results were found at 0.96 (CI 0.94-0.98). Compared to SPT and NP results ALLERG-O-LIQ shows a sensitivity of 0.88 (CI 0.83-0.92) and specificity of 0.90 (CI 0.84-0.95).




Figure 4
ROC analysis for ALLERG-O-LIQ vs. in-vivo results for seven different allergens (d1, d2, g6, t3, t4, e1 and w6) with n=287 results.

Literature

1. Kleine-Tebbe J, Breuer K, Lepp U, Vieths S, Worm M, Kunkel G, Wahn U, Lau S: **Allergen-specific IgE-values to inhalants and food allergens: Dr. Fooke ALLERG-O-LIQ versus Pharmacia CAP-System.** Allergologie 2004, 27:129-135.
2. Renz H, Becker WM, Bufe A, Kleine-Tebbe J, Raulf-Heimsoth M, Saloga J, Werfel T, Worm M: **In-vitro-Allergie-diagnostik - Positionspapier der Deutschen Gesellschaft für Allergologie und klinische Immunologie (DGAI).** Allergo J 2002, 11:492-506.
3. Mahler M, Lorenz C, Fooke M, Kleine-Tebbe J: **Comparison between skin prick test and reversed enzyme allergo sorbent test.** 20th World allergy congress, Bangkok, Thailand 2007.
4. Sun B, Mahler M, Mei J, Jing L, Zhong N: **Allergen-specific IgE to inhalant and food allergens and total IgE values in China - Comparison between two comercial immunoassays: Dr. Fooke ALLERG-O-LIQ versus Phadia ImmunoCAP® System.** J Invest All Clin Immunol 2008, 18: 319-321.
5. Reicke B, Offermann N, Fooke M: **Evaluation of a scanner based allergy lateral flow assay system for the determination of specific IgE within 20 minutes.** Oral Poster Presentation: EAACI 2017, Helsinki, Finland.

DR FOOKE