

We, the undersigned **OXY'PHARM** company, certify that our disinfection system **NOCOTECH**, as for its proven complete spectrum of activity (bactericidal, fungicidal, sporicidal and virucidal) is also efficient on **Covid-19**. **Covid-19** is part of large coronaviruses with an RNA envelope. It is similar as SARS, very contagious and highly mutant.

The precautionary methods are the same as in 2002 or 2009 : avoid hand-held or sneezing contact and disinfect surfaces as regularly as possible. Coronavirus is sensitive to H2O2, much more reactive to disinfectants than noroviruses, adenoviruses or polioviruses on which **NOCOTECH** concept gives very good results. However, it is much more contagious.

Our recommendations are as follows :
Because the precautionary principle must always prevail, we recommend daily preventive treatments of all public places. Thus, risk of contamination by surfaces will be reduced to its minimum.

ADVISED TREATMENT

Because the precautionary principle must always prevail, we recommend **daily preventive treatments** of all public places. Thus, risk of contamination by surfaces will be reduced to its minimum.

► **FREQUENCY OF USE** : Daily treatments

► **EQUIPMENT** : NOCOSPRAY® or NOCOMAX® (for big volumes)

► **PRODUCT** : NOCOLYSE® or NOCOLYS FOOD® (depending on specific needs)

► **DOSE** : 1ml/m³. Program the appliance on the volume of the room to be treated

► **DURATION OF TREATMENT** : Diffusion + 1 hour
Before re-use of the room, allow 1 hour after diffusion, then use gas detector device as per standard protocol prior to reuse of room

If a case of Covid-19 is detected, all surfaces potentially infected must be decontaminated. For that, you can make a « curative » treatment in all gathering places. We recommend to follow-up with daily 'preventive' treatments.

► **EQUIPMENT** : NOCOSPRAY® or NOCOMAX® (for big volumes)

► **PRODUCT** : NOCOLYSE ONE SHOT®

► **DOSE** : 3ml/m³ which is « volume x3 ». For the treatment of a 50m³ room program the appliance on 150m³

► **DURATION OF TREATMENT** : Diffusion + 2 hours
Before re-use of the room allow 2 hours after diffusion, then use gas detector device as per standard protocol prior to reuse of room

ROOM TREATMENT GUIDE

| ROOM VOLUME (m ³) | 10 | 20 | 30 | 40 | 50 | 100 | 150 | 200 |
|--|-------|-------|-------|-------|-------|-------|-------|-------|
| DISTRIBUTION DURATION (minute : seconde) | 03:00 | 06:00 | 09:00 | 12:00 | 15:00 | 30:00 | 45:00 | 60:00 |
| CONSUMPTION (ml) | 50 | 100 | 150 | 200 | 250 | 500 | 750 | 1000 |

Use biocidal products with care. Always read the label and product instructions before use.

HOW TO USE NOCOSPRAY


1 Connect the power cable of the computer
Set the switch to I.




2 Slide your finger on the light part of the keyboard to select the volume to be processed.
In standard mode, the selected volume default is 10m³. The selected volume range is from 10m³ to 1000m³.



3 The alleviation is progressive: if two LEDs are lit up, the intermediate value is set up.
For example, if LEDs 10 and 20 are lit up, the volume selected is 15 m³.



4 To start, press down the button in the center of the keyboard (long press).



5 A beep sounds for 15 seconds before the machine starts operating.
Leave the room before the end of 15 seconds. At the end of the broadcast, a beep sounds. The machine stops automatically.

