

# BELUNA

Negative Pressure Simulator for Breast Augmentation & Breast Firming





**CGBIO**  
CELL & GROWTH FACTOR  
BIOTECHNOLOGY

## [Negative Pressure Simulator for Breast Augmentation & Breast Firming]

Beluna puts Negative Pressure on Breast to;


- 1) Make **fat** near breast **gathered**
- 2) Improve the **flow of blood and lymph**
- 3) Help the growth of new blood vessels and cell stimulation for **regeneration of collagen and ligament**

leading to Breast Augmentation and Skin Tightening Improvement



Breast Lymphatic System






Based on the principle “that growing tissue changes its form and volume according to external stimuli”,

**Beluna increases the space by expanding skin/fat tissue using negative pressure.**




## Negative Pressure Wound Therapy

### \*Mechanism



- Gather Bloodstream > Induce Microangiopoiesis  
> Provide nutrition path for Cell proliferation and differentiation
- Provide environment for Tissue Regeneration



Creams

Y Tube +  
Filter

Dome

Adaptor

Pump

# “BEAUTY-FULL CHANGE”

Like as Crescent to Full Moon



## Better Effect

- **Powerful Motor But Compact & Silent**  
5.5L/min  
25~55mmHG
- **Wear Quickly (3sec)**
- **3 Different Mode**  
(Continuous / Cyclic / Massage)



## Ergonomic Design

- **3D Human Body**  
Curvy design
- **Stable wear fit with Female curvy design**
- **Wear less painfully**



## 1:1 Management


- **Professional consultation**
- **App for Beluna**
- **Motivate to wear Beluna everyday**



## Cream for Beluna


- **Moisture Lotion type**
- **Contain Volufiline**
- **Maintain Moisture and Increase Skin Tightening before/after Beluna wear**

## Continuous Mode




1. For Engraftment Rate Increase after Fat Graft.

## Cyclic Mode




1. For Temporary Breast Volume Increase
2. Prior to Breast Augmentation

## Massage Mode



1. After Breast Reconstruction
2. After Breast Augmentation
3. For Breast Skin Tightening Improvement



Gather the Tissue around Breast

☞ Stable and Less painful Wear

☞ Entire Breast Stimulation





- ✓ Wireless - Bluetooth connect
- ✓ Check real time Wear Condition with AI Guide
- ✓ Consultation without time and place limit
- ✓ Quick Troubleshooting

THANK YOU

Make Today, Change Tomorrow

**CGBIO**  
CELL & GROWTH FACTOR  
BIOTECHNOLOGY