

# QUIKRESPONSE®

# DIGITAL OVULATION TEST

## ONE-STEP DIGITAL HOME OVULATION PREDICTOR

WOMEN'S HEALTH



- ✓ Predicts ovulation 24 - 36 hours in advance
- ✓ Fast results in 3 minutes
- ✓ Decisive yes + / no - readout
- ✓ Available in 5, 7 and 10 kit configurations
- ✓ One-Step Test Procedure

**\*OVER 99% Accuracy at ovulation hormones**

“With Syntron’s new easy-to-read Digital Ovulation Test there’s no need to interpret confusing test results. The digital LCD display tells you in words whether you’re ovulating or not.”

### ORDERING INFORMATION:

Product Description	Packaging Option
D.O.V. - MultiStick 5 day Test Kit	10480-5
D.O.V. - MultiStick 7 day Test Kit	10480-7
D.O.V. - MultiStick 10 day Test Kit	10480-10

cGMP/ISO 13485  
CE CERTIFIED

# DIGITAL HOME OVULATION TEST

**OneStep Rapid Home Test Device for the Qualitative Detection of Human Luteinizing Hormone (LH)**

## URINE TEST

### CLEAR RESULTS


**YES +** LH Surge

**NO -** No LH Surge

- **Easy to use** – Just 1 simple step
- **Easy to read** – Clear digital results displayed

### Choose The Test Method You Prefer

Simply hold the absorbent tip directly in your urine stream for 5 seconds. Or collect your urine in a clean cup and dip the absorbent tip.



Test Stick with Wide Collection Area for Simple Sampling | Digital Reading Device | Effortless Test Stick Ejector Button | Clear Digital Display

Convenient Thumb Grip for Easy Handling

For in vitro diagnostic use. Not to be taken internally. Store in a dry place below 30 degrees C or 86 degrees F. Do not freeze. Keep out of reach of children. Not for contraceptive use or gender selection. Read instructions before use.

Certain medical conditions may affect the reliability of this test for predicting ovulation. For details see the enclosed instruction insert.

Read instructions before use. Over 99% accurate in detecting the LH surge in laboratory studies. (Based on directions for use.)

### The Digital Ovulation Predictor is easy to use and easy to read

- This Rapid, easy to use test gives you reliable results.
- Simply insert the Test Stick into Digital Reader and hold the tip of the test in your urine stream for 5 seconds.
- Clear digital results will be displayed.
- If the digital displays "YES+", you should ovulate within 24-36 hours.
- Predicts ovulation 24-36 hours in advance for early, confident family planning.
- The Digital Ovulation Predictor is a highly accurate home test that helps you predict ovulation, the days when you are most able to become pregnant.
- In clinical studies, this test has proven over 99% accurate in detecting the increase in Luteinizing Hormone ("LH surge")
- This surge usually occurs 24-36 hours before ovulation, so you can plan ahead to maximize your chances for conception.

**You may directly urinate onto the wick of the device or dip the wick of the device into a cup of urine. Within a few minutes a "YES+" or "NO-" or "?" will appear on the display screen.**

**NO-** : No LH Surge

**YES+** : LH Surge

**?** : Invalid



**SYNTRON**  
Bioresearch, Inc.

2774 Loker Ave. West | Carlsbad, CA 92010

Tel: (760) 930-2200 | Fax: (760) 930-2212 | E-mail: [marketing@syntron.net](mailto:marketing@syntron.net)

For a list of all Syntron Bioresearch, Inc. products, visit our website: [www.syntron.net](http://www.syntron.net)

Contact Syntron Bioresearch, Inc. for more information about our Women's Health Products. Call, fax, or e-mail us your order.

© 2009 Syntron Bioresearch, Inc. All rights reserved.