



## FULL FACE CPAP MASK



### ELBOW CONNECTOR

- **QUICK RELEASE**  
Quick release design allows patients to disconnect and connect to the tubing with minimum effort, increasing patient mobility.
- **MONITORING**  
Two standard 4.3mm ports are located on the connector to fulfill the demands of patient monitoring.
- **DOUBLE ROTATION**  
360 degrees rotatable connection between mask and connector, and between connector and tubing minimizes interruptions to patients.
- **EXHALATION VENT**  
Well positioned exhalation vent allows CO<sub>2</sub> to be released easily with less effort for patients.



### FRAME & MASK

- **OPTIMAL FIT**  
The reinforcing frame incorporates with the silicone mask to prevent deformation. The soft-touch cushion helps seal the mask without irritating the face.
- **EASY CLEAN**  
The connection between elbow connector and frame can be disassembled easily to allow effective cleansing.
- **SAFETY**  
Safe vent design ensures the safety of patients.
- **STERILIZATION**  
High quality medical grade materials are autoclavable to 121 °C (remove headgear before autoclave).



### HEADGEAR

- **SIMPLE CLIP**  
The clip-in design allows patients to lock the headgear into the frame and remove it efficiently.
- **GOOD FIT**  
The adjustable strap made with flexible fabric is adaptable to various patients. The permeable mesh on the backside of headgear increases patient comfort.

# NON-INVASIVE MASK



## ELBOW CONNECTOR

- **QUICK RELEASE**  
Quick release design allows patients to disconnect and connect to the tubing with minimum effort, increasing patient mobility.
- **MONITORING**  
Standard 4.3mm port is located on the connector to fulfill the demands of patient monitoring.

## FRAME & MASK

- **OPTIMAL FIT**  
The reinforcing frame incorporates with the silicone mask to prevent deformation. The soft-touch cushion helps seal the mask without irritating the face.
- **EASY CLEAN**  
The connection between elbow connector and frame can be disassembled easily to allow effective cleansing.
- **STERILIZATION**  
High quality medical grade materials are autoclavable to 121 °C (remove headgear before autoclave).

## HEADGEAR

- **SIMPLE CLIP**  
The clip-in design allows patients to lock the headgear into the frame and remove it efficiently.
- **GOOD FIT**  
The adjustable strap made with flexible fabric is adaptable to various patients. The permeable mesh on the backside of headgear increases patient comfort.