



HIGH-DEFINITION CAMERA DRAPE

The next level of camera drapes.

- Suitable for high-definition camera systems without compromising the HD image quality.
- Closed-end offers camera full protection from fluids and acts as a reliable sterile barrier.
- Allows quick and easy scope changes during procedures without compromising sterility.
- Broad-edged ring ensures ease of use and prevents contamination of the sterile field.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
ECN150-250CAP	Size 150mm x 2.5m, Ring aperture 85mm	50 units



TRAINING
VIDEOS
AVAILABLE

01 Attach the endoscope
Simply push the endoscope into the drape holder.

02 Attach the camera
Connect the camera onto the camera drape holder and extend the drape over the cable by pulling back the ring. Any excess drape is retained inside the ring.

QUICK CHANGE CAMERA DRAPE

- Simple and secure 2-step fixing method.
- Allows quick and easy scope changes during procedures without compromising sterility.
- Closed-end offers cameras full protection from fluids and acts as a reliable sterile barrier.
- Broad-edged ring reduces risk of contaminating the sterile field.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
ECE150-250C	Size 150mm x 2.5m, Ring aperture 85mm	50 units

TRAINING
VIDEOS
AVAILABLE



01 Attach the camera
Simple push-fit. The drape can then be extended over camera and cable by simply pulling back the ring. Any excess drape is retained inside the ring.

02 Attach the endoscope
Simple and secure snap-on fitting.

**TRAINING
VIDEOS
AVAILABLE**



Product code	Dimensions	Unit of sale
ECN150-250	with two adhesive tape closures Size 150mm x 2.5m, Ring aperture 85mm	50 units
ECN150-250S	elastic collar with two twist-lock spring clips and one adhesive tape closure Size 150mm x 2.5m, Ring aperture 85mm	50 units
ECN150-250TK	elastomeric end section with two adhesive tape closures Size 150mm x 2.5m, Ring aperture 85mm	50 units

CONCENTRIC RING STYLE

- Ring style drapes ensure easy and efficient application without compromising sterility.
- Drape is fully contained within the ring and can be simply pulled back over camera and cable.
- Any excess material remains securely inside the ring.
- Available with 3 different closures.
- Sterile, single-use.



TELESCOPIC STYLE

- Folded cuff to facilitate application without compromising sterility.
- Sterile, single-use.

Product code	Dimensions	Unit of sale
EOD130-250T	with perforated tip and two adhesive tape closures Size 130mm x 2.5m	25 units
EOD165-250TK	with elastomeric end section and one adhesive tape closure Size 165mm x 2.5m, Aperture 10mm dia	25 units
EOD130-250TB	with cardboard sleeve Size 130mm x 2.5m	25 units

