



MEMBER OF THE NYNOMIC GROUP

Concepts

Balloon capper

The original cap of a sample tube or vial is not always available for automatic recapping, or a secondary tube must be closed for archiving. In these cases, the most cost-effective and quick automatable decapping option is required. Many approaches are based on pre-produced disposable plastic closures or a closure with an aluminium foil. Our approach combines the advantages of both solutions, is very cost-effective and offers a very secure closure of the sample container.

The idea is to produce a hollow stopper directly in the device, which is placed in the cap like a compressed balloon and then released to tightly close the tube.

m-u-t GmbH
Am Marienhof 2
22880 Wedel
Germany

☎ +49 4103 9308 - 0
☎ +49 4103 9308 - 99
✉ info@mut-group.com
🌐 www.mut-group.com