




MTBDR-RIF 9G Membrane Kit™




Global distribution of estimated TB incidence by rate and absolute number, 2014. The size of each bubble is proportional to the size of the country's population. High burden countries are shown in red.



Source: 2015 Global Tuberculosis Report (WHO)



TB Screening + Best Drug Selection



MTBDR-RIF 9G Membrane Kit™



After PCR, simultaneous detection of MTB and Rifampin-Resistance are completed in just 30min!

Manual System

PCR



Loading & Washing






Analysis



Automation System

BMT 9G-1000



Analytical sensitivity

MTB	Limit of Detection
rpoB	10 copies/test