


Hygienic Protection ECG Belts

ECG Belts are now available with the integrated Sanitized® antimicrobial hygiene function for optimal protection.

How does it work?


The Sanitized® antimicrobial silver protects the ECG Belt.

The silver additive is a scientifically undisputed substance, whose effect on bacteria is well known. The Sanitized® technology uses a glass ceramic capsule. This ensures that the antibacterial protective treatment is firmly anchored to the ECG Belt and thus achieves long-lasting protection.



The silver ion locates the bacteria via the moisture and deactivates it.

Silver ions cannot diffuse from the ECG Belt thanks to the glass ceramic Sanitized® technology.



This is the effect silver has on bacteria:

1. Cell membrane is destabilized
2. Respiration is prohibited
3. Food (nutrient) intake is impeded
4. Cell division is inhibited

sanitized.com
Protection

