

Bio Neuro Sensor

Electroencephalogram (EEG) Sensor

Brain
Monitoring

BIS Compatible



Specifications

BME – 2002 / Adult Quattro sensor



Patient Weight.....Adult > 30kg
 Application Site.....Forehead
 Duration of UseMaximum 24 hours
 Latex Content.....Does not contain natural rubber latex
 Shelf Life.....1 year
 Operational Temperature.....-10℃~+40℃
 Storage Temperature10℃~+40℃
 Humidity.....0%~85%
 Atmosphere.....86kPa-106kPa

BME - 2003 / Pediatric sensor



Patient Weight Pediatric >20kg
 Application Site.....Forehead
 Duration of UseMaximum 24 hours
 Latex Content.....Does not contain natural rubber latex
 Shelf Life.....1 year
 Operational Temperature.....-10℃~+40℃
 Storage Temperature10℃~+40℃
 Humidity.....0%~85%
 Atmosphere.....86kPa-106kPa



Bio Medical Technologies Co.,Ltd.

101 Megu Bldg, 21-10 dongho-ro 15-gil, Jung-gu, Seoul, KOREA

Phone: 82-2-2234-0780, Fax: 82-2-2234-0748 E-mail: bmtech@bm-tech.co.kr



SC-02 (REV.00)