


# VEGAPLAX<sup>®</sup>

Advanced Plasma Technology for Medical Purpose



## Specialized for the Treatment of

- \* Plasma Arc for Blepharoplasty
- \* Plasma Lifting by Delivering Plasma into Skin Tissue

## Application Treatment Depending Power Level

- **Level 1**  
Blepharoplasty, Ring-block Anesthesia, Skin Spots, Small Soft Warts, Active Acne, Wrinkles, Periumbilical Stretch Marks, Xanthelasma, Small Neoformations or Tattoos
- **Level 2**  
Active Acne, Wrinkles, Small Keloids, Small Ruby Angiomas, Face Lifting, Significant Xanthelasma with Thickness
- **Level 3**  
Large Moles, Keloids, Medium / Large Ruby Angiomas, Seborrhic Warts, Large Fibromas

## Specifications

- **Output Frequency** : 70kHz
- **Power Level** : 3 Steps (0.7W, 1.0W, 2.0W)
- **Output Mode** : Continuous & Pulse
- **Pulse Mode Frequency** : 10, 20, 40, 100, 200Hz
- **Life Time** : Fully charged, Max 3 hours
- **Power** : Internal Li-Polymer Rechargeable Battery 7.4V, 1,000mAh
- **This product is guaranteed for 1 year.**
- **General life of recharging battery is 500 times of charging and discharging.**



*Acmedi* KOREA

# Before & After

## Blepharoplasty



24 hours later



72 hours later



4 days later



5 days later



7 days later




Right after Treatment

## Wrinkle Improvement



Before



After



Before



After

## Scar Removal



Before



After

## Acne & Blemish Treatment



Before



After