

Maccura Rapid test for COVID-19

Severe Acute Respiratory Syndrome Coronavirus 2(SARS-CoV-2) IgM/IgG Antibody Assay Kit by Colloidal Gold Method



Product feature

IgM and IgG 2-in-1

2 results from 1 test, which can simultaneously determine negative/positive status of both IgM and IgG.

The guarantee of material optimization

The optimized protein fragments N and S of SARS-CoV-2 can improve detection capacity.

Quick interpretation of results


Under normal circumstances, the results can be observed within 10-15 minutes. After 15 minutes, the results are invalid.

Strong positive results can be observed in 150 seconds.


Simple and convenient operation

Evaluating results without the aid of instruments

Operation procedure



Result interpretation



Product information

Product name	Specification	Catalog No.
Severe Acute Respiratory Syndrome Coronavirus 2(SARS-CoV-2) IgM/IgG Antibody Assay Kit by Colloidal Gold Method	20 Tests	IM4102079