


CATHCHECK™

Instant epidural catheter verification with CompuWave™ displaying the pulsatile waveform found in the epidural space.

A new, modern technique for checking your catheter placement in real time.



*Catheter sold separately

Benefits

EPIDURAL CATHETER VERIFICATION SYSTEM

- » Check catheters in seconds, not minutes
- » Saves the anesthesia provider time, enabling them to make clinical decisions faster regarding ineffective epidural catheters and the next course of treatment
- » The anesthesia provider may flush with sterile anesthetic or sterile saline
- » Increases the anesthesia provider's daily effectiveness
- » Provides confirmation to the anesthesia provider post-catheter placement
- » Reduces the amount of time the patient is in pain, due to an ineffective catheter

	Catheter Dosing	CompuFlo® CathCheck™
Connects directly to the patient's catheter	X	X
Flush with anesthetic	X	X
Flush with sterile saline		X
Time required for patient response	40 minutes	2 minutes
Recognizes the presence of pulsatile waveform when the epidural space is accessed with CompuWave™		X