



Product Construction Sheet

FLX068885 - FLEXcon® dermaFLEX™ TT 200 BLACK H-502 200 POLY H-9 44PP-8

GENERAL DESCRIPTION:

This is a 2.0 mil (51 micron) permanent, pressure-sensitive, acrylic adhesive with organic polar salt and a conductive matrix supported on one side with a tight-release, clear, smooth, polyester release liner (200 POLY H-9), and on the other side with a free-release, poly-coated, one-side, natural-kraft-paper, release liner (44PP-8).

The adhesive has been tested for biocompatibility (skin irritation, sensitization and cytotoxicity) by an independent testing agency specializing in biocompatibility referencing ISO-10993. The adhesive also meets the requirements of AAMI/ANSI EC-12*

PRODUCT DATA	VALUE	TEST METHOD
PHYSICAL PROPERTIES		
Thickness, mils (µm)	200 POLY H-9 Liner: 2.0 (51) +/- 10% Adhesive: 2.0-2.2 (51-56) +/- 0.1 (3) 44PP-8 Liner: 3.0 (76) +/- 10 %	ASTM D 3652 (Modified for use with non-tape product)
ADHESION PROPERTIES**		
	Average	
180° Peel Stainless Steel	13 oz/in (143 N/m)	ASTM D 903 15 min dwell at RT
Expected Shear	2 hours	ASTM D 3654 Method A (1 hr. dwell, 1 sq. in, 4 lb. load)
Tack	340 g	ASTM D 2979
ELECTRICAL PROPERTIES**		
Z-Direction Resistance	<20 (Ohms/in ²)	ASTM D 257
Surface Resistivity	<100 (Ohms/□)	ASTM D 257
STORAGE STABILITY	Two years stored at 70°F (21°C) and 50% RH in the original packaging.	

*DISCLAIMER: FLEXcon MAKES NO WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO ANY IMPLIED WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, COURSE OF PERFORMANCE, OR TRADE USAGE. This product is a device component intended solely for further processing, manufacturing, incorporation into a finished device, and is further intended only for use as a device component on intact skin (i.e., not intended for use in the presence of open sores or wounds nor for applications inside the body). Customer is solely responsible for determining whether this component product is fit for Customer's intended use, including, without limitation, incorporation into Customer's finished device(s), and for performing any additional testing that may be necessary to support Customer's intended use. Customer is solely responsible for ensuring that its use of the component product is in compliance with the applicable laws and for obtaining any necessary clearances or approvals from its use of the component product and any finished device that incorporates the component product. Full biocompatibility results are available upon request.

Product Performance and Suitability

All of the descriptive information, the typical performance data, and recommendations for the use of FLEXcon products shall be used only as a guide and do not reflect the specification or specification range for any particular property of the product. Furnishing such information is merely an attempt to assist you after you have indicated your contemplated use and shall in no event constitute a warranty of any kind by FLEXcon. All purchasers of FLEXcon products shall be responsible for independently determining the suitability of the material for the purpose for which it is purchased. No distributor, salesman, or representative of FLEXcon is authorized to give any warranty, guaranty, or make any representation in addition or contrary to the above.

** For guidance only, not to be used for setting specifications. The data presented are based on three or less production runs. Product construction sheets are updated as necessary, without notice. Please contact FLEXcon to request a copy of most recent version.

dermaFLEX™ is a trademark of FLEXcon Company, Inc.
