

## Product Datasheet

### SARS-CoV-2 Pool

### Human Plasma

#### Application:

Pool of Human Plasma, SARS-CoV-2 IgG, is used for production purposes like production of standards, calibrators or controls for IVDs or EQAS.

<b>Product No.:</b>	DSPO 4.1.9.1.1
<b>Volume</b>	1.0 l individual volumes on request
<b>Source Material</b>	true human CPD plasma, pooled
<b>Parameter</b>	SARS-CoV-2 IgG (IgM values on availability)
<b>Details given</b>	titer, volume, matrix, test kit, test manufacturer
<b>Bottling</b>	100.0 ml
<b>Test Samples</b>	available
<b>Disease State Testing</b>	Each individual donation has been tested with methods cleared in compliance with European Directive 98/79/EC and found negative for HBsAg, anti-HIV I+II, anti-HCV, HIV I RNA, HCV RNA, HBV DNA, HCV RNA, Syphilis.
<b>Materials Safety</b>	All material of human origin is to be regarded as potentially infectious. Therefore standard precautions are necessary, if handling the material.  All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals.
<b>Storage</b>	≤ -15°C
<b>Minimum Shelf Life</b>	48 months
<b>Manufacturing Quality System</b>	in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485.
<b>Pricing</b>	Please inquire for your individual quote.

# Contact & Ordering Information



## Central Contact

**+49 (0)3302 55 199 0**  
**info@inventdiagnostica.de**

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

## Contact for your Order

**sales@inventdiagnostica.de**

