

EasyHip Cementless Cup

Wide range of sizes
 • 42 - 60 mm, 62 mm on request

Triple locked inserts
 Secure fixation of polyethylene insert to shell by a "snap-lock mechanism", conical coupling and antirotation tabs

Variable bone screw options
 • Shells with pre-drilled holes for additional screw fixation
 • Sector Hole Shell screw holes are pre-closed
 • +/-15° screw angulation possible²

Color coding for streamlined workflow



Streamlined surgical experience through simple, ergonomic instruments and color coding

Lipped Liner: 5mm raised "antilux" insert.

50/36 – Outside small, inside big
 36 mm prosthesis heads starting from 50 mm shell size for ceramic and polyethylene inserts

Latest material selection
 • Biolox® Delta Ceramic
 • Vit-E polyethylene



TPS+CaP double coating promotes bony ongrowth^{1,3}
 ➔ A highly porous 300 µm thick Titanium Plasma Spray (TPS) coating and an electrochemically deposited 15 µm thick Calcium Phosphate (CaP) layer to ensure optimal primary and secondary implant stability.

