













# MICROMETHOD

From patient to sampling, LABELIANS identifies the solution



Patient group	Sampling applications	Puncture site	Specifications Diameter Depth	Recommended product/Ref.
 Adults and children over 1 year old	Blood sampling	Fingers: • middle finger • ring finger	Blade 18G 1.8mm	 <b>LSPRESS18G1.8</b>
				 <b>LSLITE18G18</b>
	Blood glucose	Fingers: • middle finger • ring finger	Needle 30G 1,8mm	 <b>LSPRESS30G1.8</b>
 Children under 1 year old	Guthrie's test, blood sampling	Heel	Blade 0,85mm	 <b>LSBABYP11</b>
	Blood glucose	Big toe	Needle 30G 1,8mm	 <b>LSPRESS30G1.8</b>
Term child and over 3kg	Guthrie's test, blood sampling	Heel	Blade - 1mm	 <b>LSBABYP05</b>
			Blade - 2mm	 <b>LSBABYP06</b>
	Blood glucose	Big toe	Needle 30G 1,8mm	 <b>LSPRESS30G1.8</b>

\*Examples of applications, given for information purposes only.

