

PRIZMA PABPM

Technical Data

- Technology: Osilometric method
- Systolic blood pressure: 40-260 mmHg
- Diastolic blood pressure: 20-210 mmHg
- Pulse: 40-200 bpm
- Accuracy: ± 3 mmHg
- Battery: 1000 mAh Li battery, approx. 100 measurements
- Memory: up to 200 measurement results
- Dimensions: device 119 × 52 × 21 mm, cuff 26-36 cm
- Weight: 105 g (with battery)

CONTENTS OF THE KIT:

PRIZMA PABPM ambulatory blood pressure monitor,
M cuff (upper arm circumference 26-36 cm),
USB cable, and instruction manual.

Optional accessories:

S cuff (upper arm circumference 18-26 cm),
L cuff (upper arm circumference 30-43 cm)

Without tube

All in one

On The Phone

Settings and results

Clinical report

Statistics and graphics



CE0120