




VESALIUS SERIES

LaproArm-f

laparoscope manipulator

"YOUR THIRD HAND"



SUSTAINABLE SURGICAL SOLUTIONS

VESALIUS SERIES: LaproArm-f

AN "ASSISTANT" WHO DOES EXACTLY
WHAT YOU WANT AND NEVER WEAKENS!

The devices from the VESALIUS series work like the ideal assistant. The systems are built according to our unique RCM principle. They are compact, easy to use and can be operated with one hand. Due to their unique compact structure, they do not interfere with the other instruments during manipulation.

VESALIUS LaproArm-f : laparoscope manipulator

The VESALIUS LaproArm-f is a mechanical laparoscope manipulator that is operated simply and with one hand. According to the weight and the position of the laparoscope the surgeon can adjust at any time the friction braking force. If necessary, the laparoscope can always be disconnected from the device and moved freely. Re-coupling is done equally quickly. By combining a standard drape and sterilizable accessories, the device can be used frequently and cheaply on a daily basis for virtually all laparoscopic surgeries.

+ Stable images

Even in the most difficult positions, the use of the VESALIUS LaproArm-f offers a stable and fixed image.

+ Bionic arm

This device from the VESALIUS series behave like an arm with the laparoscope holder acting as the forearm. It is fixed in volume, and does not pivot. As a result, it will not unexpectedly interfere with other instruments during manipulation.

+ Sterility

The VESALIUS LaproArm-f is covered with a standard drape. The sterilizable accessories (adaptor and instrument holder) are autoclavable. As a result, there are virtually no running costs.

+ Ergonomic and efficient

The surgeon gains more control over her/his operation. Even in extreme positions, the image is retained. Highly trained assistants can perform demanding tasks more efficiently.

+ Intensive use


The durable device is installed in a matter of minutes. As a result it can provide support for virtually all laparoscopic procedures on a daily basis.

+ Wide range of applications


The VESALIUS LaproArm-f can be used in almost all types of laparoscopic surgery: general surgery, gynecology and urology.



Movement: left - right



Movement: up - down



Warning: Before using our products, please read carefully the sections on hazards, warnings, safety instructions and indications in the product manual and also follow the warning symbols on the products. **Note:** The above information is subject to change without notice. The standard description, accessories and technical specifications are in accordance with our price list, not the photographs of the machines presented in the catalog. It is forbidden to reproduce or copy, in whole or in part, texts, drawings and images without prior written permission.

**HEPHAESTUS
SURGICAL**

HEPHAESTUS Surgical
Dubbeekstraat 112
3200 Aarschot
info@hephasurgical.com

© by METAL AARSCHOT nv
first published Belgium 2018 / VESALIUS LaproArm-f / v. 0.1 - 03.2018 /
Printed in Belgium / Subject to technical changes

HEPHAESTUS Surgical™ and VESALIUS LaproArm-f™ are trademarks of
METAL AARSCHOT nv



Designed and Made in Belgium

www.hephasurgical.com