

# iM70 Patient Monitor



**iM70**  
Patient Monitor



## About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Patient Monitoring
- OB/GYN
- Ultrasound Imaging
- Point-of-Care Testing
- *In-Vitro* Diagnostics
- Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



**Global Headquarters:**  
Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen  
518122 P.R. China | +86.755.26898326 | [www.edan.com](http://www.edan.com) | [info@edan.com](mailto:info@edan.com)

**U.S. and Canada inquiries:**  
EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126  
+1.858.750.3066 | [www.edandiagnostics.com](http://www.edandiagnostics.com) | [edan-info@edandiagnostics.com](mailto:edan-info@edandiagnostics.com)

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.



ENG-PM-iM70-V1.0-20200311



# iM70

## Patient Monitor



12.1"

Maximum 13 Waveforms



No-fan Design Night Mode



Pitch Tone Thermal Recorder



Shortcut Menu Ethernet Printer




### Standard Parameters

- 3/5-lead ECG, HR, RESP, EDAN SpO<sub>2</sub>, EDAN NIBP, PR, 2-TEMP

### Optional Parameters


- 12-lead ECG, 2-IBP(with ICP support), C.O., EDAN G2 CO<sub>2</sub> (Sidestream), Respiration CO<sub>2</sub> (Mainstream and Sidestream), MASIMO Anesthesia Gas/O<sub>2</sub>, Nellcor OxiMax™ SpO<sub>2</sub>, SunTech® NIBP



### Reliable Algorithms

- iSEAP™ ECG algorithm optimized for arrhythmia detection, pacemaker detection, and HR measurement
- SEMIP® 12-lead diagnosis algorithm with 208 kinds of analysis results
- iMAT™ SpO<sub>2</sub> algorithm with outstanding motion resistance and low perfusion resistance performance
- iCUFFS™ NIBP algorithm optimized for cardiac patients, hypertensive patients, and neonatal patients

### Multiple Display Modes



### End tidal Carbon Dioxide for Intubated/Non-Intubated Patients



#### EDAN G2 CO<sub>2</sub> (Sidestream)

- Superior water trap design for accurate monitoring
- iCARB™ algorithm with intelligent CO<sub>2</sub> pseudo wave identification technology
- Multiple sampling accessories as options for adult, child and neonate patients

#### Respiration CO<sub>2</sub> (Mainstream/Sidestream)

- Plug & play module design
- Dehumidification tube instead of water trap
- Low sampling rate of 50ml/min suitable for all types of patients

### MASIMO Anesthetic Gas/O<sub>2</sub> (Mainstream/Sidestream)




- Unique Plug & Play module for mainstream anesthetic gas monitoring
- The no moisture sampling line design suitable for all clinical applications
- A stable low flow of 50ml/min for all types of patients



Nurse Call



VGA Output



Defib Proof & Sync



ESU Proof



Bed Rail Hook



Wall Mount



Rolling Stand

### Built-in Non-volatile Memory:

A single piece of patient data contains:

1200 NIBP Measurement Review	120 h Trend Review	200 Alarm Review	48 h Full disclosure Wave
---------------------------------	-----------------------	---------------------	------------------------------