



deepskin

IONIC SKIN CARE SYSTEM

The whitening and wrinkle management technology, widely known as Vital Ion (iontophoresis) in dermatology and beauty salons, has been miniaturized.

Please experience the home aesthetic system with which people can get the effects of skin whitening and wrinkle improvement anywhere because it is compatible with various mask packs and it administers a functional serum deep into the skin.



Main Features



Effective

- It makes a functional essence directly penetrate into the skin in an activated state
- The skin penetration rate increases up to 8.6 times at 15 minutes after using it (preclinical results of the device performance)



Safe

- DeepSkin causes little irritation to the skin and there are no issues of addiction or side effects
- DeepSkin was verified through clinical trials and there is a patent for it
- South Korea KC, Europe CE, US UL, FCC Certified



Economical

- DeepSkin can be used conveniently at home
- DeepSkin is economical through saving time and cost (Less time, less cost)




Convenient

- Automatic management without touching it through connecting the mask to the main body
- The design is as small as a makeup compact
- Convenient portability owing to light weight



Iontophoresis Effect





Iontophoresis

As one of the transdermal drug delivery systems, this technology involves injecting ionized nutrients deep into the dermal layer of the skin through the repulsive power and attractive force of electric charges using micro-currents.

Skin penetration amount increases
8.6 times

DeepSkin Preclinical results of the device performance
<Sookmyung Women's University Pharmaceutical Research Institute>

Clinical Results



(Crow's feet area, Derma TOP)

※ Results of skin clinical tests of the device through DeepSkin mask packs <every day for 12 weeks, 25 women>

※ There may be differences depending on individual skin conditions or the use environment.

Hicare Ionic Infusion Mask

Whitening/wrinkle improvement dual functional product



✓ Skin elasticity improves with Collagen Essence

- Contain hydrolyzed collagen derived from fish
- Supply moisture and nutrition to the skin with a molecular structure similar to the skin
- Maintain the firm skin and improve the skin tone through reinforcing skin elasticity

✓ Use Berry Complex Ingredients


- Contain 5 types of berry extracts; Strawberry, Blueberry, Raspberry, Blackberry and Acai Berry
- Skin Soothing, Sebum Control and Skin Moisturizing Effects
- Contain anthocyanin to prevent skin aging

✓ Highly-concentrated 100% Tencel Cupra Mask

- An eco-friendly material which enhances the strength of Tencel to the soft adhesion of Cupra
- Improve moisture retention by delivering the essence with excellent adhesion (Essence absorption power more than 10 times compared with the weight of the sheet)



Main Components



Specification

Product Name (Model Name)	DeepSkin (HE-100S)
Item	Skincare device
Applicable Technology	Iontophoresis Input Technology (Pulse type)
Frequency	110Hz
Operation Time	Maximum 30 Minutes (5-Minute Unit)
Time of Use	Max. 30 min. a day recommended (Initial setting: 15 min.)
Intensity Control	1~10 Stages (Initial setting: 10 stages)
Size	Breadth(65mm) X Length(65mm) X Height(25mm)
Weight	44g(batteries excluded)
Power	4.5V(AAA 1.5V X 3)