




OKM8800

Central Monitoring System (CMS)



OKM8800 Central Monitoring System

Minimum hardware configuration

Operation system	Windows XP, Windows 10 (32 bit), and Windows 10 (64 bit)
Display	17-inch, 19-inch, 21-inch color CRT or LCD monitor
Resolution	1280 x 1024 pixels or above, with aspect ratio of 4:3, 5:4, 16:10
Network card	PC with integrated 10M/100M adaptive Ethernet NIC
Graphics card	256M or above independent video memory
CPU	2.8GHz or above Intel dual-core CPU
Memory	2G or above
Hard disk	500G or above
Sound card	16 bits

Review

Dynamic short trend	2 hours
Trend review	Trend table and trend graph, mass storage
Waveform review	Mass storage
NIBP review	Mass storage
Alarm review	Mass storage of Alarm events including the all parameter values at the alarm moment and waveform of 16s before and after the alarm moment.
Historical review	Up to 20,000 historical patient data.

Calculation

Drug calculation	Up to 750 calculations for a single patient
Hemodynamics calculation	Up to 750 calculations for a single patient
Oxygenation calculation	Up to 750 calculations for a single patient
Ventilation calculation	Up to 750 calculations for a single patient
Renal calculation	Up to 750 calculations for a single patient

Printing

Print paper/format	A4, A5, B5 (195*245)mm, B5 (182*257)mm, B5 (176*250)mm, C5, 16K
Print contents	Real-time waveform, 7-lead ECG, 12-lead ECG, waveform review, fetal waveform review, trend table, trend graph, alarm event list, alarm event review, patient information.

Data import and export

File format	Binary file
Storage medium	Mobile hard disk, flash drive, floppy disk
Contents	Patient information, waveform data, parameter trend data, alarm events, drug calculations, hemodynamics calculations, oxygenation calculations, ventilation calculations, renal calculations.

OKUMAN designs and produces medical devices for OR, NICU, and ICU. The product range covers anesthesia devices, patient monitoring, infant incubators, ventilators, medical pendants, operation lights and tables. All devices designed with specific patient data managing system to link-in to multiple medical database.

Specifications subject to change without notice

Okuman Medikal Sistemler A.Ş.

Add: Kazım Karabekir Cad. 95/95 06060 İskitler Ankara/Turkey

Tel: +90 312 384 05 20

Web: www.okuman.com.tr

Fax: +90 312 384 19 75

E-mail: info@okuman.com.tr