

Rigid Endoscopy Trolley SMOVE

ENDOMED
S Y S T E M S



The **SMOVE** trolley is designed for rigid endoscopy systems and for the specific devices.

It is easy to move and build with a modular system, which makes adding or removing of related accessories and devices a lot easier. The color of the device is a special medical grade.

- Solid & Stable
- Highly mobile
- Modular system
- Medical Paint
- Clip break feature
- Colored elements
- Easy clean design

SMOVE

EndoMed Systems Trolley

Item No.	Description
648000	EndoMed Trolley 'SMOVE'
648001	Shelf 450x 360 mm
648002	Shelf 450 x 360 mm with handle
648003	Shelf with key board tray and mousepad
648004	Drawer 450 x 360 mm
648005	Push handle
648006	Endoscope holder
648007	Twin monitor holder fixed
648008	Twin monitor holder adjustable
648009	Isolating transformer 600VA
648010	Isolating transformer 1200VA

EndoMed Systems GmbH
88214 Ravensburg, Germany
www.endomed.com

Tel.: +49 (0)751 35978 0
Fax: +49 (0)751 35978 102
email@endomed.com



Rigid Endoscopy Trolley SMOVE



Technical Specifications

Device	SMOVE
Beam Height	1172 mm
Width and Depth	450 x 360 mm
Paint	Light grey
Wheels	4 breakable Twin wheels Ø 100mm
Material	High strength Aluminium
Loading Capacity	80 Kg
Wheel points	610 x 620 mm
Total Height	1415 mm
Weight	46 Kg
Delivery	
Standard Configuration (648000)	<ul style="list-style-type: none"> • 1x Basic frame with 4 lockable twin castors Ø 100 m • 2x Shelves and one drawer • 1x Swivel and tilt monitor holder, adjustable up to 14 Kg • 1x Endoscope holder with clamp • 1x 6 equipotential bonding connectors • 1x On /off switch • 2x Rear cable winders
Options	<ul style="list-style-type: none"> • Keyboard drawer with mousepad (left/right) • Push handle • Endoscope holder • Twin monitor holder fixed • Twin monitor holder adjustable • Isolating transformer 600VA • Isolating transformer 1200VA
Environmental Conditions	
Temperature (Operation)	+10°C to +40°C
Temperature (Transport and Storage)	-20° C to +45°C
Max. relative humidity	During operation: 70% During storage and transport: 70%

EndoMed Systems GmbH
88214 Ravensburg, Germany
www.endomed.com

Tel.: +49 (0)751 35978 0
Fax: +49 (0)751 35978 102
email@endomed.com

